



Intensity SL



* Intensity Toric



* Available in
preloaded configuration

	Intensity SL	Intensity Toric *
Platform	C-Loop	Plate Haptic
Overall Diameter	13 mm	11 mm (>16D) 11.5 mm (≤16D)
Placement	Capsular Bag	
Optic Diameter	6 mm	
Power Range	10 to 30D (0.5D increments)	
Addition	Add powers: +3 and +1.5 in IOL plane for near and intermediate; and two additional focus intensifying diffractive orders providing the lens a greater efficacy	
Cylinder Range for Powers 10.0 to 20.0	-	1.0, 1.5, 2.25, 3.0
Cylinder Range for Powers 20.5 to 30.0	-	1.0, 1.5, 2.25, 3.0, 3.75, 4.5
Optic Design	Pentafocal Posterior surface: Aspheric - Diffractive Anterior surface: Spherical/Toric Pupil aperture optimized	
360° Continuous Square Edge	Yes	
Haptic Angulation	5°	0°
Material	Hydrophilic acrylic	
Refractive Index	1.46 (@ 35°C)	
SRK/T A-Constant Non-Contact Biometry	118.4	117.45
SRK/T A-Constant Contact Biometry	118.06	117.11
Spherical Aberration	-0.13μ	
Light Filtration	Natural Yellow Violet Filter	